RapidVit™ Omni

Vitrolife

US

en: Indication for use

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EUROPE

en: Indication for use

Media for vitrification of oocytes through to blastocyst stage embryos.

Contraindications

RapidVit™ Omni contains gentamicin. Do not use in patients with known hypersensitivity/allergy to the component.

Product Description

SUPPLEMENTED WITH HSA

RapidVit Omni contains three solutions for the vitrification of oocytes through to blastocyst stage embryos. The solutions consist of a MOPS-buffered medium containing gentamicin as an antibacterial agent, and human serum albumin.

Vitri 1™ Omni contains no cryoprotectants.

Vitri 2™ Omni contains ethylene glycol and propanediol as cryoprotectants.

Vitri 3™ Omni contains ethylene glycol and propanediol as cryoprotectants and sucrose.

For use after warming to +37 °C in ambient atmosphere.

The summary of safety and clinical performance can be found at www.vitrolife.com

Storage instructions and stability

Store dark at +2 to +8 °C.

RapidVit Omni is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles can be used for up to two weeks after first opening, use aseptic technique and minimize the time outside the refrigerator. Record opening date on the bottle. Discard excess media no later than two weeks after first opening.

US: Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

Directions for use

The product shall be used by an IVF professional.

The patient target group is an adult or reproductive-age population that undergoes IVF treatment or fertility preservation.

The following is the general procedure for using RapidVit Omni.

Note: Timing with vitrification procedures is critical, ensure you follow the protocol precisely.

Vitrification

Place 1 ml of each of the following media into separate wells of a multi-well plate and warm to 37 °C in ambient atmosphere.

• Vitri 1 Omni

• Vitri 2 Omni

• Vitri 3 Omni

All manipulations of the oocytes and embryos are carried out at 37 °C (on the heated stage) in ambient atmosphere. It is very important to keep the temperature of all vitrification solutions at 37°C at all times. Deviations from 37°C will alter the permeability of cryoprotectants, which may compromise oocyte and embryo survival.

NOTE: The recommended volumes should not be changed. Failure to use the correct volume of media may result in osmolality changes, which could cause suboptimal oocyte and embryo survival.

Oocytes:
Following egg retrieval oocytes should be held in an incubator in equilibrated G-1™ PLUS at 37°C at 6% CO_2. Approximately two hours after retrieval denude the oocytes. If oocytes are not vitrified immediately after denudation, they should be placed in G-1 PLUS until the time of vitrification.

Transfer the oocytes into Vitri 1 Omni. Oocytes should remain in the solution for at least 5 minutes, but a maximum of 20 minutes.

Move an appropriate number of oocytes into Vitri 2 Omni. The oocytes will be exposed to this solution for 2-5 minutes. The oocytes will tend to float to the surface, if so, collect and replace them to the bottom of the dish.

The oocytes are expected to re-expand to their original volume within 5 minutes. If re-expansion is not observed within 5 minutes this could indicate that the temperature of the medium is not at 37°C. Temperature will affect the rate at which the cryoprotectants enter the oocyte and therefore the observed re-expansion.

Prepare the cryodevice for use.

When the oocytes have fully re-expanded make two 50 µl droplets of Vitri 3 Omni on a non-toxic surface, preferably a culture dish. Making small droplets of Vitri 3 Omni enables easy loading onto the cryodevice.

Embryos:
Blastocysts may be collapsed prior to dehydration procedure using an appropriate biopsy pipette. Puncture of the blastocyst should be done at a location furthest away from the inner cell mass.

Transfer the embryos into Vitri 1 Omni. Embryos should remain in the solution for at least 5 minutes, but a maximum of 20 minutes.

Move an appropriate number of embryos into Vitri 2 Omni. The embryos will tend to float to the surface, if so, collect and replace them to the bottom of the dish.

The embryos will be exposed to this solution for 2 minutes.

Towards the end of the two minutes, make two 50 µl droplets of Vitri 3 Omni on a non-toxic surface, preferably a culture dish. Making small droplets of Vitri 3 Omni enables easy loading onto the cryodevice.

Use of storage device

Use a legally marketed storage device indicated for use in human oocyte and embryo vitrification procedures, in order to assure sufficient cooling and warming rates.

Use a closed storage system to prevent the potential risk of viral or other contamination of samples.

Perform the actual vitrification and warming according to the *Instructions for Use* for the storage device.

NOTE: The 50 µl droplets can only be used once.

Transfer the oocytes/embryos in a minimal volume of Vitri 2 Omni to avoid dilution of the droplet.

Transfer the oocytes/embryos into the first Vitri 3 Omni droplet. Empty the pipette outside the droplet and prime the pipette again from the second droplet. Immediately transfer the oocytes/embryos to the second Vitri 3 Omni droplet.

NOTE: The total time from transferring the oocytes/embryos into the droplet of Vitri 3 Omni until the oocytes/embryos are vitrified should be between 30-45 seconds.

Immediately vitrify the oocytes/embryos according to the cryodevice instruction.

Continue storage according to laboratory practice.

Specifications

Aseptically filtered

Mouse Embryo Assay, 1-cell [% embryos developed to expanded blastocyst at 96 hours]	≥ 80
Bacterial endotoxins (LAL assay) [EU/mL]	< 0.5
pH tested	
Osmolality tested	

LOT-specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions

The long-term safety of vitrification and/or blastocyst collapse on children born following this method of oocyte or embryo cryopreservation procedure has not been established.

The risk of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife´s IVF media, have not been determined and are uncertain.

Discard product if bottle integrity is compromised. Do not use RapidVit Omni if it appears cloudy.

RapidVit Omni contains human serum albumin.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for antibodies to HIV, HBc, HCV, and HTLV I/II and non-reactive for HbsAg, HCV RNA and HIV-1 RNA and syphilis. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Not for injection.

Discard the product according to standard clinical practice for medical hazardous waste when the procedure is finished.

Caution: Federal (US) law restricts this device to sale by or on the order of a physician or practitioner trained in its use (Rx only).

Symbols glossary

* Medical devices – Symbols to be used with information to be supplied by the manufacturer (ISO 15223-1:2021).

Symbol	Title	Reference and symbol ID
	Consult instructions for use	*Ref no: 5.4.3
	Caution	*Ref no: 5.4.4
	Manufacturer	*Ref no: 5.1.1
	Date of manufacturer	*Ref no: 5.1.3
	Use-by date	*Ref no: 5.1.4
LOT	Batch code	*Ref no: 5.1.5
REF	Catalogue number	*Ref no: 5.1.6
STERILE A	Sterilized using aseptic processing techniques	*Ref no: 5.2.2
	Single sterile barrier system	*Ref no: 5.2.11
+2°C +8°C	Temperature limit	*Ref no: 5.3.7
	Keep away from sunlight (light sources)	*Ref no: 5.3.2
	Contains a medicinal substance	*Ref no: 5.4.7
	Contains human blood or plasma derivatives	*Ref no: 5.4.6
MD	Medical Device	*Ref no: 5.7.7
CE 2460	CE Mark (Conformité Européenne)	European conformity (CE) mark with Notified Body identification number as defined in Regulation (EU) 2017/745 Annex V

Caution: Federal (US) law restricts this device to sale by or on the order of a physician or practitioner trained in its use (Rx only).

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