RapidWarm™ Omni

Vitrolife

US

en: Indication for use

Media for warming of vitrified oocytes through to blastocyst stage embryos.

Contraindications

RapidWarm™ Omni contains gentamicin. Do not use in patients with known hypersensitivity/allergy to the component.

Product Description

SUPPLEMENTED WITH HSA

RapidWarm Omni contains four solutions for the warming of vitrified oocytes through to blastocyst stage embryos. The solutions consist of a MOPS- buffered medium containing gentamicin as an antibacterial agent, and human serum albumin.

Warm 1™ Omni contains sucrose

Warm 2™ Omni contains sucrose

Warm 3™ Omni contains sucrose

Warm 4™ Omni contains no sucrose

For use after warming to +37 °C in ambient atmosphere.

The summary of safety and clinical performance can be found at www.vitrolife.com

Storage instructions and stability

Store dark at +2 to +8 °C.

RapidWarm™ Omni is stable until the expiry date shown on the container and the LOT-specific Certificate of Analysis.

Media bottles can be used for up to two weeks after first opening, use aseptic technique and minimize the time outside the refrigerator. Record opening date on the bottle. Discard excess media no later than two weeks after first opening.

US: Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

Directions for use

The product shall be used by an IVF professional.

The patient target group is an adult or reproductive-age population that undergoes IVF treatment or fertility preservation.

The following is the general procedure for using RapidWarm Omni.

Note: Timing with warming is critical, ensure you follow the protocol precisely.

Warming:

Place 1 ml of each of the following media into separate wells of a multi-well plate and warm to 37 °C in ambient atmosphere:

- Warm 1 Omni
- Warm 2 Omni
- Warm 3 Omni
- Warm 4 Omni

NOTE: The recommended volumes should not be changed. Volume changes will affect temperature control in the first warming solution as well as osmolality, which may result in suboptimal oocyte or embryo survival.

It is very important to keep the temperature of all warming solutions at 37°C at all times. If more than one cryodevice is to be warmed in the same dish, make sure that the temperature reaches 37°C after each warming procedure.

All manipulations of the oocytes or embryos are carried out at 37°C (on the heated stage) in ambient atmosphere.

Remove the cryodevice containing the vitrified oocytes or embryos from the cryostorage container. Follow the warming instructions for the specific cryodevice used.

Immediately, transfer oocytes or embryos into Warm 1 Omni. It is important to use 1ml of Warm 1 Omni to ensure the temperature of the medium is not perturbed by successive cryodevice submersions.

Allow the oocytes or embryos to fall from the device and sink to the bottom. Leave for 1 minute.

Transfer the oocytes or embryos into Warm 2 Omni and let the oocytes or embryos remain in the solution for 3 minutes.

Transfer the oocytes or embryos into Warm 3 Omni and let the oocytes or embryos remain in the solution for 5 minutes.

Transfer the oocytes or embryos into Warm 4 Omni and let the oocytes or embryos remain in the solution for 5 to 10 minutes.

Following warming, oocytes or embryos should be held in an incubator in equilibrated G-1™ PLUS at 37°C and 6% CO2 for minimum 2 hours before continuing culture according to standard laboratory practice.

Storage instructions and stability

Store dark at +2 to +8 °C.

RapidWarm™ Omni is stable until the expiry date shown on the container and the LOT-specific Certificate of Analysis.

Media bottles can be used for up to two weeks after first opening, use aseptic technique and minimize the time outside the refrigerator. Record opening date on the bottle. Discard excess media no later than two weeks after first opening.

US: Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

Directions for use

The product shall be used by an IVF professional.

The patient target group is an adult or reproductive-age population that undergoes IVF treatment or fertility preservation.

The following is the general procedure for using RapidWarm Omni.

Note: Timing with warming is critical, ensure you follow the protocol precisely.

Warming:

Place 1 ml of each of the following media into separate wells of a multi-well plate and warm to 37 °C in ambient atmosphere:

- Warm 1 Omni
- Warm 2 Omni
- Warm 3 Omni
- Warm 4 Omni

NOTE: The recommended volumes should not be changed. Volume changes will affect temperature control in the first warming solution as well as osmolality, which may result in suboptimal oocyte or embryo survival.

It is very important to keep the temperature of all warming solutions at 37°C at all times. If more than one cryodevice is to be warmed in the same dish, make sure that the temperature reaches 37°C after each warming procedure.

All manipulations of the oocytes or embryos are carried out at 37°C (on the heated stage) in ambient atmosphere.

Remove the cryodevice containing the vitrified oocytes or embryos from the cryostorage container. Follow the warming instructions for the specific cryodevice used.

Immediately, transfer oocytes or embryos into Warm 1 Omni. It is important to use 1ml of Warm 1 Omni to ensure the temperature of the medium is not perturbed by successive cryodevice submersions.

Allow the oocytes or embryos to fall from the device and sink to the bottom. Leave for 1 minute.

Transfer the oocytes or embryos into Warm 2 Omni and let the oocytes or embryos remain in the solution for 3 minutes.

Transfer the oocytes or embryos into Warm 3 Omni and let the oocytes or embryos remain in the solution for 5 minutes.

Transfer the oocytes or embryos into Warm 4 Omni and let the oocytes or embryos remain in the solution for 5 to 10 minutes.

Following warming, oocytes should be held in an incubator in equilibrated G-1™ PLUS at 37°C at 6% CO_2. The oocytes should be washed several times through this medium.

Perform ICSI 1-2 hours after warming. Continue culture according to standard laboratory practice.

Following warming, embryos should be held in an incubator in equilibrated G-1 PLUS, G-TL™,G-2™ PLUS or EmbryoGlue® at 37°C at 6% CO_2 according to standard laboratory practice.

Specifications

Aseptically filtered	
Mouse Embryo Assay, 1-cell [% embryos developed to expanded blastocyst at 96 hours]	≥ 80
Bacterial endotoxins (LAL assay) [EU/mL]	< 0.5
pH tested	
Osmolality tested	

LOT-specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions

The long-term safety of vitrification and/or blastocyst collapse on children born following this method of oocyte or embryo cryopreservation procedure has not been established.

The risk of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife´s IVF media, have not been determined and are uncertain.

Discard product if bottle integrity is compromised. Do not use RapidWarm Omni if it appears cloudy.

RapidWarm Omni contains human serum albumin.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for antibodies to HIV, HBc, HCV, and HTLV I/II and non-reactive for HbsAg, HCV RNA and HIV-1 RNA and syphilis. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Not for injection.

Discard the product according to standard clinical practice for medical hazardous waste when the procedure is finished.

Symbols glossary

* Medical devices – Symbols to be used with information to be supplied by the manufacturer (ISO 15223-1:2021).

Symbol	Title	Reference and symbol ID
	Consult instructions for use	*Ref no: 5.4.3
	Caution	*Ref no: 5.4.4
	Manufacturer	*Ref no: 5.1.1
	Date of manufacturer	*Ref no: 5.1.3
	Use-by date	*Ref no: 5.1.4
LOT	Batch code	*Ref no: 5.1.5
REF	Catalogue number	*Ref no: 5.1.6
STERILE A	Sterilized using aseptic processing techniques	*Ref no: 5.2.2
	Single sterile barrier system	*Ref no: 5.2.11
+2°C / +8°C	Temperature limit	*Ref no: 5.3.7
	Keep away from sunlight (light sources)	*Ref no: 5.3.2
	Contains a medicinal substance	*Ref no: 5.4.7
	Contains human blood or plasma derivatives	*Ref no: 5.4.6
MD	Medical Device	*Ref no: 5.7.7
CE 2460	CE Mark (Conformité Européenne)	European conformity (CE) mark with Notified Body identification number as defined in Regulation (EU) 2017/745 Annex V

Caution: Federal (US) law restricts this device to sale by or on the order of a physician or practitioner trained in its use (Rx only).

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